



Recombinant Human Calmodulin-like protein 4 (CALML4)

Product Code	CSB-MP836218HU
Abbreviation	CALML4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96GE6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAEHLLPGP PPSLADFRLE AGGKGTERGS GSSKPTGSSR GPRMAKFLSQ DQINEYKECF SLYDKQQRGK IKATDLMVAM RCLGASPTPG EVQRHLQTHG IDNGELDFS TFLTIMHMQI KQEDPKKEIL LAMLMVDKEK KGYVMASDLR SKLTSLGEKL THKEVDDLFR EADIEPNGKV KYDEFIHKIT LPGRDY
Source	Mammalian cell
Target Names	CALML4
Protein Names	Recommended name: Calmodulin-like protein 4 Alternative name(s): Serologically defined breast cancer antigen NY-BR-20
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.