



# Recombinant Human Probable ATP-dependent RNA helicase DDX27 (DDX27), partial

<b>Product Code</b>	CSB-YP836220HU
<b>Abbreviation</b>	DDX27
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96GQ7
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DNMLADLGLIGTIGEDDEVPVEPESDSGDEEEEGPIVLGRRQKALGKNRSADF NPDFVFTEKEGTYDGSWALADVMSQLKKKRAATTLDEKIEKVRKKRKTEDKEA KSGKLEKEKEAKEGSEPKEQEDLQENDEEGSEDEASET DYSSADENILTKADT LKVKDRKKKKKKGQEAGGFFEDASQYDENLSFQDMNLSRPLLKAITAMGFKQ PTPIQKACIPVGLLGKDICACAATGTGKTAAFALPVLERLIYKPRQAPVTRVLVL VPTRELGIQVHSVTRQLAQFCNITTCLAVGGGLDVKSQEAALRAAPDILIATPGRL IDHLHNCPSFHLSSIEVLILDEADRMLDEYFEEQMKEIIRMCSHHRQTMLFSAT MTDEVKDLASVSLKNPVRIFVNSNTDVAPFLRQEFIRIRPNREGDREAIVAALLT RTFTDHSVMLFTQTKKQAHRMHILLGLMGLQVGELHGNLSQTQRLEALRRFKD EQIDILVATDVAARGLDIEGVKTVINFTMPNTIKHYVHRVGR TARAGRAGRSVS LVGEDERKMLKEIVKAAKAPVKARILPQDVILKFRDKIEKMEKD VYAVLQLEAEE KEMQQSEAQINTAKRLLKKGKEAVVQEPERSWFQTKEERKKEKIAKALQEFDL ALRGKKRKKKFMKDAKKKGEMTAEERSQFEILKAQMFAERLAKRNRRAKRAR AMPEEEPVRGPAKKQKQGKKS VFDEELTNTSKKALKQYRAGPSFEERKQLGL PHQRRGGNFKSKSRYKR
<b>Source</b>	Yeast
<b>Target Names</b>	DDX27
<b>Protein Names</b>	Recommended name: Probable ATP-dependent RNA helicase DDX27 EC= 3.6.4.13 Alternative name(s): DEAD box protein 27
<b>Expression Region</b>	30-794aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.