



Recombinant Human Protein lin-37 homolog (LIN37)

Product Code	CSB-YP836222HU
Abbreviation	LIN37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96GY3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFPVKVKVEK SELEMAKARN QLDAVLQCLL EKSHMDRERL DEEAGKTPSD THNKDCSIAA TGKRPSARFP HQRRKKRREM DDGLAEGGPQ RSNTYVIKLF DRSVDLAQFS ENTPLYPICR AWMRNSPSVR ERECSPSSPL PPLPEDEEGS EVTNSKSRDV YKLPPPTPPG PPGDACRSRI PSPLQPEMQG TPDDEPSEPE PSPSTLIYRN MQRWKRIRQR WKEASHRNQL RYSESMKILR EMYERQ
Source	Yeast
Target Names	LIN37
Protein Names	Recommended name: Protein lin-37 homolog Alternative name(s): Antolefinin
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein expressed in the eye.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.