



Recombinant Human LIM domain-containing protein ajuba (AJUBA)

Product Code	CSB-BP836229HU
Abbreviation	AJUBA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96IF1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MERLGEKASR LLEKFGRRKG ESSRSGSDGT PGP GKGR LSG LGGPRKSGPR GATGGPGDEP LEPAREQGSL DAERNQRGSF EAPRYEGSFP AGPPPTRALP LPQSLPPDFR LEPTAPALSP RSSFASSAS DASKPSSPRG SLLLDGAGAG GAGGSRPCSN RTSGISMGYD QRHGSPLPAG PCLFGPPLAG APAGYSPGGV PSAYPELHAA LDRLYAQRPA GFGCQESRHS YPPALGSPGA LAGAGVGAAG PLERRGAQPG RHSV TGYGDC AVGARYQDEL TALLRLTVGT GGREAGARGE PSGIEPSGLE EPPGPFVPEA ARARMREPEA REDYFGTCIK CNKGIYQSN ACQALDSL YH TQCFVCCSCG RTL RCKAFYS VNGSVYCEED YLFSGFQ EAA EKCCVCGHLI LEKILQAMGK SYHPGCFRCI VCNKCLDGIP FTVDFSNQVY CVTDYHKNYA PKCAACGQPI LPSEGCE DIV RVISMDRDYH FECYHCEDCR MQLSDEEGCC CFPLDGHLLC HGCHMQRLNA RQPPANYI
Source	Baculovirus
Target Names	AJUBA
Protein Names	Recommended name: LIM domain-containing protein ajuba
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.