



Recombinant Human TNFAIP3-interacting protein 3 (TNIP3)

Product Code	CSB-YP836243HU
Abbreviation	TNIP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96KP6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAHFVQGTSR MIAAESSTEH KECAEPSTRK NLMNSLEQKI RCLEKQRKEL LEVNQQWDQQ FRSMKELYER KVAELKTKLD AAERFLSTRE KDPHQQRKRD DRQREDDRRQ DLTRDRLQRE EKEKERLNEE LHELKEENKL LKGKNTLANK EKEHYECEIK RLNKALQDAL NIKCSFSEDC LRKSRVEFCH EEMRTEMEVL KQQVQIYEED FKKERSDRER LNQEKEELQQ INETSQSQLN RLNSQIKACQ MEKEKLEKQL KQMYCPPCNC GLVFHLQDPW VPTGPGAVQK QREHPPDYQW YALDQLPPDV QHKANGLSSV KKVHP
Source	Yeast
Target Names	TNIP3
Protein Names	Recommended name: TNFAIP3-interacting protein 3 Alternative name(s): A20-binding inhibitor of NF-kappa-B activation 3 Short name= ABIN-3 Listeria-induced gene protein
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.