



Recombinant Human Protein FAM167A (FAM167A)

Product Code	CSB-EP836244HU
Abbreviation	FAM167A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96KS9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSVPQIHVEE VGAEEGAGAA APPDDHLRSL KALTEKLRLE TRRPSYLEWQ ARLEEHTWPF PRPAAEPQAS LEEGERGGQE PLLPLREAGQ HPPSARSASQ GARPLSTGKL EGFQSIDEAI AWLRKELTEM RLQDQQLARQ LMRLRGDINK LKIEHTCRLH RRMLNDATYE LEERDELADL FCDSPLASSF SLSTPLKLI GTKMNINSRR FSLC
Source	E.coli
Target Names	FAM167A
Protein Names	Recommended name: Protein FAM167A
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.