



Recombinant Human Protein FAM81B (FAM81B)

Product Code	CSB-EP836252HU-B
Abbreviation	FAM81B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LP2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQLQFLGTLA SSEKRKKSQR LFFKNIKSTK NKAGKASIMS SDTNVNKSAS PTATAEEQPV EPDGPLPGSD NNQEKVRLS PAKMSTKNST DLVEYVDKSH AFLPIIPNTQ RGQLEDRLNN QARTIAFLLQ QAFRIKEDIS ACLQGTHGFR KEESLARKLL ESHIQTITSI VKKLSQNI EI LEDQIRARDQ AATGTNFAVH EINIKHLQGV GDLRGRVARC DSSIVKLSGD IHLFRQEHRQ IEKAIQEFVP ALETLSKNLD MKVMQLLGKI ETASPEQTSN LKMVQGDYRH EMNLLEFKFH SLSSNLYEEV ENNKKWTENQ FLKYRKDHLG HINECLKVLQ EKLEKSENKM EEKLLQLSSK VENFINTQKQ ETQLSKVKHM ENKLSKKMEQ MEKQIWGELE TMQNEYQSGF KSIHDSLSSL QQIQKTKMDL EKYKVQKDLK KLQRKIVELQ EV
Source	E.coli
Target Names	FAM81B
Protein Names	Recommended name: Protein FAM81B
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.