



Recombinant Human Protein SPATIAL (SPATIAL)

Product Code	CSB-MP836256HU
Abbreviation	SPATIAL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M53
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MATDVQLADY PLMSPKAELK LEKKSGRKPR SPRDSGPQKE LVIPGIVDFE RIRRALRTPK PQTGTYCFG RLSHHSFFSR HHPHPQHVTI IQDLTGKPVC VVRDFPAPLP ESTVFSGCQM GIPTISVPIG DPQSNRNPQL SSEAWKKELK ELASRVAFLT KEDELKKKEK EQKEEPLREQ GAKYSAETGR LIPASTRAVG RRRSHQGQQS QSSSRHEGVQ AFLLDQDQELL VLELLCRILE TDLLSAIQFW LLYAPPKEKD LALGLLQTAV AQLLPQPLVS IPTEKLLSQL PEVHEPPQEK QEPPCSQSPK KTKISPFTKS EKPEYIGEAQ VLQMHSSQNT EKKTSKPRAE S
Source	Mammalian cell
Target Names	TBATA
Protein Names	Recommended name: Protein SPATIAL Alternative name(s): Stromal protein associated with thymii and lymph node homolog
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.