



Recombinant Human UPF0536 protein C12orf66 (C12orf66)

Product Code	CSB-EP836259HU-B
Abbreviation	C12orf66
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MD2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGESIPLAAP VPVEQAVLET FFSHLGIFSY DKAKDNVEKE REANKSAGGS WLSLLAALAH LAAAEKVYHS LTYLGQQLGG QSFFSRKDSI RTIYTSLHNE LKKVVTGRGA LGGTAPHVEE LLSHLSEQLC FFVQARMEMA DFYEKMYTLS TQKFINAEEL VGLLDAIMKK YSSRFHHPIL SPLESSFQLE VDVLCHELLKA QAQVSEWKFL PSLVNLHSAH TKLQTWGQIF EKQRETKKHL FGGQSQKAVQ PPHLFLWLMK LKNMLLAKFS FYFHEALSRQ TTASEMKTTLT AKANPDFFGK ISSFIRKYDA ANVSLIFDNR GSESFQGHGY HHPHSYREAP KGVDQYPAVV SLPSDRPVMH WPNVIMIMTD RTSDLNSLEK VVHFYDDKVQ STYFLTRPEP HFTIVIIFES KKSERDSHF I SFLNEVSLAL KNPKV FASLK PGAKG
Source	E.coli
Target Names	C12orf66
Protein Names	Recommended name: UPF0536 protein C12orf66
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.