



# Recombinant Human UPF0536 protein C12orf66 (C12orf66)

<b>Product Code</b>	CSB-MP836259HU
<b>Abbreviation</b>	C12orf66
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96MD2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGESIPLAAP VPVEQAVLET FFSHLGIFSY DKAKDNVEKE REANKSAGGS WLSLLAALAH LAAAEKVYHS LTYLGQKLGQ QSFFSRKDSI RTIYTSLHNE LKKVVTGRGA LGGTAPHVEE LLSHLSEQLC FFVQARMEMA DFYEKMYTLS TQKFIAEEL VGLLDAIMKK YSSRFHHPIL SPLESSFQLE VDVLCHELLKA QAQVSEWKFL PSLVNLHSAH TKLQTWGQIF EKQRETKKHL FGGQSQKAVQ PPHLFLWLMK LKNMLLAKFS FYFHEALSRQ TTASEMKTTLT AKANPDFFGK ISSFIRKYDA ANVSLIFDNR GSESFQGHGY HHPHSYREAP KGVDQYPAVV SLPSDRPVMH WPNVIMIMTD RTSDLNSLEK VVHFYDDKVQ STYFLTRPEP HFTIVIIFES KKSERDSHF I SFLNEVSLAL KNPKV FASLK PGAKG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	C12orf66
<b>Protein Names</b>	Recommended name: UPF0536 protein C12orf66
<b>Expression Region</b>	1-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.