



Recombinant Human Noelin-3 (OLFM3)

Product Code	CSB-BP836274HU
Abbreviation	OLFM3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96PB7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	QTLPSLV GLNTRLSTP DTLTQISPKE GWQVYSSAQD PDGRRICTVV APEQNLCSR D AKSRQLRQLL EKVQNMSQSI EVLNLRTQRD FQYVLKMETQ MKGLKAKFRQ IEDDRKTLMT KHQELKEKM DELLPLIPVL EQYKTDAKLI TQFKEEIRNL SAVLTGIQEE IGAYDYEELH QRVLSLETRL RDCMKKLTGC KLMKITGPVT VKTSGTRFGA WMTDPLASEK NNRVWYMSY TNNKIVREYK SIADFVSGAE SRTYNLPFKW AGTNHVYNG SLYFNKYQSN IIKYSFDMG RVLAQRSLEY AGFHNVPYPT WGGFSDIDLM ADEIGLWAVY ATNQNAGNIV ISQLNQDTLE VMKSWSTGYP KRSAGESFMI CGTLYVTNSH LTGAKVYYSY STKTSTYEYD DIPFHNQYFH ISMLDYNARD RALYAWNNGH QVLFNVTLFH IIKTEDDT
Source	Baculovirus
Target Names	OLFM3
Protein Names	Recommended name: Noelin-3 Alternative name(s): Olfactomedin-3 Optimedin
Expression Region	24-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.