



# Recombinant Human CCA tRNA nucleotidyltransferase 1, mitochondrial (TRNT1)

<b>Product Code</b>	CSB-MP836282HU
<b>Abbreviation</b>	TRNT1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96Q11
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TEGLKSLTE LFKENHELR IAGGAVRDLL NGVKPQDIDF ATTATPTQMK EMFQSAGIRM INNRGEKHGT ITARLHEENF EITTLRIDVT TDGRHAEVEF TTDWQKDAER RDLTINSMFL GFDGTLFDYF NGYEDLKNKK VRFVGHAKQR IQEDYLRILR YFRFYGRIVD KPGDHDPEL EAI AENAKGL AGISGERIWW ELKKILVGNH VNHLIHLIYD LDVAPYIGLP ANASLEEFDK VSKNVDGFSP KPVTLASLF KVQDDVTKLD LRLKIAKEEK NLGLFIVKNR KDLIKATDSS DPLKPYQDFI IDSREPDATT RVCELLKYQG EHCLLKEMQQ WSIPFPVSG HDIRKVGISS GKEIGALLQQ LREQWKKSGY QMEKDELLEY IKKT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TRNT1
<b>Protein Names</b>	Recommended name: CCA tRNA nucleotidyltransferase 1, mitochondrial EC=2.7.7.72 Alternative name(s): Mitochondrial tRNA nucleotidyl transferase, CCA-adding mt CCA-adding enzyme mt tRNA CCA-diphosphorylase mt tRNA CCA-pyrophospho
<b>Expression Region</b>	42-434
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.