



Recombinant Human Acidic repeat-containing protein (ACRC)

Product Code	CSB-MP836285HU
Abbreviation	ACRC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QF7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDGCKKELPR LQEPEEDED C YILNVQSSSD DTSGSSVARR APKRQASCIL NVQSRSGDTS GSSVARRAPK RQASSVVID SDSDEECHTH EEKAKLLEI NSDDESPECC HVKPAIQEPP IVISDDDNDD DNGNDLEVPD DNSDDSEAPD DNSDDSEAPD DNSDDSEAPD DNSDDSEAPD DNSDDSDVPD DNSDDSSDDN SDDSSDDNSD DSDVPDDKSD DSDVPDDSSD DSDVPDDSSD DSEAPDDSSD DSEAPDDSSD DSEAPDDSSD DSEAPDDSSD DSEASDDSSD DSEASDDSSD DSEAPDDKSD DSDVPEDKSD DSDVPDDNSD DLEVPVPAED LCNEGQIASD EEELVEAAAA VSQHDSSDDA GEQDLGENLS KPPSDPEANP EVSERKLPTE EEPAPVVEQS GKRKSKTKTI VEPPRKRQTK TKNIVEPPRK RQTKKNIVE PLRKRKAKTK NVSVTPGHKK RGPSKKKPGA AKVEKRKTRT PKCKVPGCFL QDLEKSKKYS GKNLKRNKDE LVQRIYDLFN RSVCDKKLPE KLRIGWNNKM VKTAGLCSTG EMWYPKWRRF AKIQIGLKVC DSADRIRDTL IHEMCHAASW LIDGIHDSHG DAWKYYARKS NRIHPELPRV TRCHNYKINY KVHYECTGCK TRIGCYTKSL DTSRFICAKC KGS LVMVPLT QKDGTRIVPH V
Source	Mammalian cell
Target Names	GCNA
Protein Names	Recommended name: Acidic repeat-containing protein
Expression Region	1-691
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.