



Recombinant Human Protein FRG2-like-1 (FRG2B)

Product Code	CSB-EP836287HU-B
Abbreviation	FRG2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QU4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGKGNEDPDL HCSSIQCSTD QPPFQQISFT EKGSDKKPF KEKGKTAFSH SSEKHIQRQG SEPNNPKENS EETKLKAGNS TAGSEPESSS YRENCRKRM SSKDSCQDTA GNCPEKECSL SLNKKSRST PVHNSEIQT CDAHHRGRSR ACTGRSKRHR SRALGVQTPS IRKSLVTSVR AMSEAVYQDL AQVWAQQIHS PLTCEQLTLL TRLRGPLCAQ VQTLYSMATQ AAYVFPAESW LVPATLPGPG ESALDREAHP FPGQEITETV SGSDEAKL
Source	E.coli
Target Names	FRG2B
Protein Names	Recommended name: Protein FRG2-like-1 Alternative name(s): FSHD region gene 2 protein family member B HSA10-FRG2
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.