



# Recombinant Human Nucleus accumbens-associated protein 1 (NACC1)

<b>Product Code</b>	CSB-YP836291HU
<b>Abbreviation</b>	NACC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96RE7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAQTLQMEIP NFGNSILECL NEQRLQGLYC DVSVVVKGHA FKAHRAVLAA SSSYFRDLFN NSRSVVVLP AAVQPQSFQQ ILSFCYTGRL SMNVGDQFLL MYTAGFLQIQ EIMEKGTEFF LKVSSPSCDS QGLHAEAEAPS SEPQSPVAQT SGWPACSTPL PLVSRVKTEQ QESDSVQCMP VAKRLWDSGQ KEAGGGGNGS RKMAKFSTPD LAANRPHQPP PPQQAPVVAA AQPAVAAGAG QPAGGVAAAG GVVSGPSTSE RTSPGTSSAY TSDSPGSYHN EEDEEEDGGE EGMDEQYRQI CNMYTMYSMM NVGQTAEKVE ALPEQVAPES RNRIRVRQDL ASLPAELINQ IGNRCHPKLY DEGDPSEKLE LVTGTNVYIT RAQLMNCHVS AGTRHKVLLR RLLASFFDRN TLANSCGTGI RSSTNDPRRK PLDSRVLHAV KYQCQNFAPN FKESEMNAIA ADMCTNARRV VRKSWMPKVK VLKAEDDAYT TFISETGKIE PDMMGVEHGF ETASHEGEAG PSAEALQ
<b>Source</b>	Yeast
<b>Target Names</b>	NACC1
<b>Protein Names</b>	Recommended name: Nucleus accumbens-associated protein 1 Short name= NAC-1 Alternative name(s): BTB/POZ domain-containing protein 14B
<b>Expression Region</b>	1-527
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.