



Recombinant Human Tribbles homolog 1 (TRIB1)

Product Code	CSB-EP836296HU-B
Abbreviation	TRIB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96RU8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRVGPVRSAM SGASQPRGPA LLFPATRGVP AKRLLDADDA AAVAAKCPRL SECSSPPDYI SPPGSPCSPQ PPPAAPGAGG GSGSAPGPSR IADYLLLPLA EREHVSRALC IHTGRELRCK VFPIKHYQDK IRPYIQLPSH SNITGIVEVI LGETKAYVFF EKDFGDMHSY VRSRKRLREE EAARLFKQIV SAVAHCHQSA IVLGDCLKRK FVFSTEERTQ LRLESLEDTH IMKGEDDALS DKHGCPAYVS PEILNTTGTY SGKAADVWSL GVMLYTLLVG RYPFHSDSDPS ALFSKIRRGQ FCIPEHISPK ARCLIRSLLR REPSERLTAP EILLHPWFES VLEPGYIDSE IGTSDQIVPE YQEDSDISSF FC
Source	E.coli
Target Names	TRIB1
Protein Names	Recommended name: Tribbles homolog 1 Short name= TRB-1 Alternative name(s): G-protein-coupled receptor-induced gene 2 protein Short name= GIG-2 SKIP1
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.