




Recombinant *Aspergillus kawachii* Probable endo-beta-1,4-glucanase D (eglD)

Product Code	CSB-YP836318APO
Relevance	Has endoglucanase activity on substrates containing beta-1,4 glycosidic bonds, like in carboxymethylcellulose (CMC), hydroxyethylcellulose (HEC) and beta-glucan. Involved in the degradation of complex natural cellulosic substrates
Abbreviation	Recombinant <i>Aspergillus kawachii</i> eglD protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96WQ9
Product Type	Recombinant Protein
Immunogen Species	<i>Aspergillus kawachii</i> (strain NBRC 4308) (White koji mold) (<i>Aspergillus awamori</i> var. <i>kawachi</i>)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	HTTVQAVWINGEDQGLGNTDDGYIRSPSPNSPVTDVTSTDMTCNVNGDQAASKTLSVKAGDVVTFEWHHSRSDSDSDIIASSHKGPVQVYMAPTAKGSNGNNWVKIAEDGYHKSSDEWATDILIANKGKHNITVPDVPAGNYLFRPEIIALHEGNREGGAQFYMECVQFKVTS DGSNELPSGVSPGVYTATDPGILFDIYNSFDSYPIPGPDVWDGSSSGSSSSGSSSSAAVSSAAAAATTSAVAATTPATQAAVEVSSSSAAAATTEAAAPVVSSAAPVQQATSAVTSQAQAAPTTFATSSKKSSKTACKNKTKSNSQVAAATSSV VAPAATSSVVPVVSASASASAGGVAKQYERCGGINHTGPTTCESGSVCKKWNPPYYQCVASQ
Research Area	others
Source	Yeast
Target Names	eglD
Protein Names	Carboxymethylcellulase D Cellulase 61A Cellulase D cel61A
Expression Region	21-408aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	41.7 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.