



Recombinant Arabidopsis thaliana E3 ubiquitin-protein ligase RING1-like (At3g19950)

Product Code	CSB-EP836568DOA-B
Abbreviation	At3g19950
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LPN7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSSGVNSTGS AAAAPEVDKM FFCYQCNQTV TISISSADP FCPICNQGFL EEYEDPNPNQ SLNFPNNSD SFFPMADPFS TLLPLIFGSS AAAPSGMDFM SLFGPSMQPQ ARSTQQNPQS DAFDPFTFLQ NHLQTLRSSG THFEFVIENH PSDPGNRMPG NFGDYFFGPG LEQLIQQLAE NDPNRYGTPP ASKSAIDALP TVKVTKDMLK SEMNQCAVCM DEFEDGSDVK QMPCKHVFHQ DCLLPWLELH NSCPVCRFEL PTDDPDYENR SQGSQGSQGDG QGSVEGQQTTP RFSIQLPWPFP RRQDGSGSGS GAPGTGGGNL ETRGEDLD
Source	E.coli
Target Names	At3g19950
Protein Names	Recommended name: E3 ubiquitin-protein ligase RING1-like EC= 6.3.2.- Alternative name(s): RING finger protein 1
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.