



Recombinant *Oryza sativa* subsp. *japonica* Auxin-responsive protein IAA6 (IAA6)

Product Code	CSB-MP836571OFG
Abbreviation	IAA6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LQ74
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MEEGSNKREG LPPQLLDLIP DEKEWKLREA LGLGRSRNAG FDGEEDKKLD LKLGLPGFIE DDEAETLRDY RLQQESPSLS LSFFPKHSKT TSSTTTTGA KRGFIDTVED KTEGYNDQKQ QARAGCGKEL AVEEMIAAVS ERKKGCCPPP PPPHGAPATP ARNRPQTQGR GAAAPVVGWP PIRSFRRNLA SSSSSKHSPE PQNDNANAKV TLTCKKNPLV KINMDGIPIG RKIDLAAYNS YDGLSSAVKQ LFHGFLQAQK DQTNAQIAQQ GADDKIFYQL LDGSGEYTLV YEDSEGDRML VGDVPWKV FV STAKRLRVLR SSELSTLIG ATARV
Source	Mammalian cell
Target Names	IAA6
Protein Names	Recommended name: Auxin-responsive protein IAA6 Alternative name(s): Indoleacetic acid-induced protein 6
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.