



Recombinant Human Putative ribosome-binding factor A, mitochondrial (RBFA)

Product Code	CSB-MP836636HU
Abbreviation	RBFA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N0V3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	KNWLKKFAS KTKKKVWYES PSLGSHSTYK PSKLEFLMRS TSKKTRKEDH ARLRALNGLL YKALTDLLCT PEVSQELYDL NVELSKVSLT PDFSACRAYW KTTLSAEQNA HMEAVLQRSA AHMRHLLMSQ QTLRNVPIV FVQDKGNAAL AELDQLLAVA DFGPRDERDN FVQNDFRDPD APQPCGTTEP TTSSSLCGID HEALNKQIME YKRRKDKGLG GLVWQQQVAE LTTQMKKGRK RAKPRLEQDS SLKSYLSGEE VEDDLDLVGA PEYECYAPDT EELEAERGGG RTEDGHSCGA SRE
Source	Mammalian cell
Target Names	RBFA
Protein Names	Recommended name: Putative ribosome-binding factor A, mitochondrial
Expression Region	42-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.