



# Recombinant Human Rho GTPase-activating protein 24 (ARHGAP24)

<b>Product Code</b>	CSB-MP836648HU
<b>Abbreviation</b>	ARHGAP24
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N264
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEENNDSTEN PQQGQGRQNA IKCGWLRKQG GFVKTWHTRW FVLKGDQLYY FKDEDETKPL GTIFLPGNKV SEHPCNEENP GKFLFEVVP GDRDRMTANH ESYLLMASTQ NDMEDWVKS IRRVIWGPFGG GIFGQKLED VRYEKRYG NR LAPMLVEQCV DFIRQRGLKE EGLFRLPGQA NLVKELQDA DCGEKPSFDS NTDVHTVASL LKLYLRELPE PVIPYAKYED FLSCAKLLSK EEEAGVKELA KQVKSLP VVN YNLLKYICRF LDEVQSYSGV NKMSVQNLAT VFGPNILRPK VEDPLTIMEG TVVVQQLMSV MISKHDCLFP KDAELQSKPQ DGVSNNNEIQ KKATMGQLQN KENNNTKDSP SRQCSWDKSE SPQRSSMNG SPTALSGSKT NSPKNSVHKL DVSRSPPLMV KKNPAFNKGS GIVTNGSFSS SNAEGLEKTQ TTPNGSLQAR RSSSLKVS GT KMGTHSVQNG TVRMGILNSD TLGNPTNVRN MSWLPNGYVT LRDNKQKEQA GELGQHNRLS TYDNVHQQFS MMNLDDKQSI DSATWSTSSC EISLPENSNS CRSSTTTCPE QDFFGGNFED PVLDPGPPQDD LSHPRDYESK SDHRVSGGRS SRATSSSDNS ETFVGNSSSN HSAHSLVSS LKQEMTKQKI EYESRIKSLE QRNLTLET EM MSLHDELDQE RKKFTMIEIK MRNAERAKED AEKRNDMLQK EMEQQFFSTFG ELTVEPRRTE RGNTIWIQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ARHGAP24
<b>Protein Names</b>	Recommended name: Rho GTPase-activating protein 24 Alternative name(s): Filamin-A-associated RhoGAP Short name= FilGAP RAC1- and CDC42-specific GTPase-activating protein of 72 kDa Short name= RC-GAP72 Rho-type GTPase-activati
<b>Expression Region</b>	1-748
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.