



Recombinant Human RIB43A-like with coiled-coils protein 1 (RIBC1)

Product Code	CSB-EP836661HU-B
Abbreviation	RIBC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N443
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MYNIKQSTDT KEAAAIEARR NREKERQNR FNVNRVVMGV DVQALNNQVG DRKRREAAER SKEAAYGTSQ VQYDVVVQML EKEEADRTRQ LAKKVQEFRE KQKQLKNGRE FSLWDPGQVW KGLPTYLSYS NTYPGPASLQ YFSGEDLDRD TRLRMQQGQF RYNLERQQQE QQQAKVDENY TDALSNQLRL AMDAQATHLA RLEESCRAAM MCAMANANKA QAAVQAGRQR CERQREQKAN LAEIHQSTS DLLTENPQVA QHPMAPYRVL PYCWKGMTPE QQAAIRKEQE VQRSKKQHR QAECTLDTEW KSQTMSSAQA VLELEEQERE LCAVFQRGLG SFNQQLANEQ KAQQDYLSV IYTNQPTAQY HQQFNTSSR
Source	E.coli
Target Names	RIBC1
Protein Names	Recommended name: RIB43A-like with coiled-coils protein 1
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.