



# Recombinant Human RIB43A-like with coiled-coils protein 1 (RIBC1)

<b>Product Code</b>	CSB-MP836661HU
<b>Abbreviation</b>	RIBC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N443
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MYNIKQSTDT KEAAAIEARR NREKERQNR FNVNRVVMGV DVQALNNQVG DRKRREAAER SKEAAYGTSQ VQYDVVVQML EKEEADRTRQ LAKKVQEFRE KQKQLKNGRE FSLWDPGQVW KGLPTYLSYS NTYPGPASLQ YFSGEDLDRD TRLRMQQGQF RYNLERQQQE QQQAKVDENY TDALSNQLRL AMDAQATHLA RLEESCRAAM MCAMANANKA QAAVQAGRQR CERQREQKAN LAEIQHQSTS DLLTENPQVA QHPMAPYRVL PYCWKGMTPE QQAAIRKEQE VQRSKKQHR QAECTLDTEW KSQTMSSAQA VLELEEQERE LCAVFQRGLG SFNQQLANEQ KAQQDYLSV IYTNQPTAQY HQQFNTSSR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RIBC1
<b>Protein Names</b>	Recommended name: RIB43A-like with coiled-coils protein 1
<b>Expression Region</b>	1-379
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.