



Recombinant Human Coiled-coil domain-containing protein 24 (CCDC24)

Product Code	CSB-EP836665HU-B
Abbreviation	CCDC24
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4L8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLRHSPSLWE LVEEHVPLRE RREVKRILGE AAVDLSLELR AEVAMLRALL QEARSSQAPS SRPISDPSSL LAPPPLLKDL LRQELRQLLQ GLRHKAIPEG RDQAQAWVQY SPRVLHFALE EPRCDLPEQE IFQMRGGGPS SGHRDLSIIK DQLNVSNIQ VARHLRGLLE EECHTLEREI LILQRCLEEE YLRPCHPSEA ALEPTLAELK EQKKAMEQEL QASVGPSCVS PNHRQRPLGS STQGLRPPLP LCGVAPLQCC LPAPPLEPYL RPRGQSATHR WGRQLQCSPR EGPASTPMSS AAPQAPA
Source	E.coli
Target Names	CCDC24
Protein Names	Recommended name: Coiled-coil domain-containing protein 24
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.