



Recombinant Human HORMA domain-containing protein 2 (HORMAD2)

Product Code	CSB-MP836683HU
Abbreviation	HORMAD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N7B1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MATAQLSHCI TIHKASKETV FPSQITNEHE SLKMOVKKLFA TSISCITYLR GLFPESYGE RHLDDLKSLKI LREDKKCPGS LHIIRWIIQGC FDALEKRYLR MAVLTLYTDP MGSEKVTEMY QKFKYTKEG ATMDFDSHSS STSFESGTNN EDIKKASVLL IRKLYILMQD LEPLPNNVVL TMKLHYNAV TPHDYQPLGF KEGVNSHFLF FDKEPINVQV GFVSTGFHSM KVKVMTEATK VIDLENNLFR ENSTTEIAHQ GLDCDEEEEC NDHIQRMNFV CSQQSSECSR KKRKVSEPVK VFIPNRK
Source	Mammalian cell
Target Names	HORMAD2
Protein Names	Recommended name: HORMA domain-containing protein 2
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.