



Recombinant Human Transcription cofactor vestigial-like protein 2 (VGLL2)

Product Code	CSB-EP836691HU
Abbreviation	VGLL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N8G2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSCLDVMYQV YGPPQPYFAA AYTPYHQKLA YYSKMQEAQE CNASPSSSGS GSSSFSSQTP ASIKEEEGSP EKERPPEAEY INSRCVLFTY FQGDIVSVVD EHFSTRALSQP SSYSPSCTSS KAPRSSGPWR DCSFPMSQRS FPASFWNSAY QAPVPPPLGS PLATAHSELP FAAADPYSPA ALHGHLHQGA TEPWHHAPH HAHPHHPYAL GGALGAQAAP YRPAAVHEV YAPHFDPRYG PLLMPAASGR PARLATAPAP APGSPPELS GKGEPAAGAAW AGPGGPFASP SGDVAQGLGL SVDSARRYSL CGASLLS
Source	E.coli
Target Names	VGLL2
Protein Names	Recommended name: Transcription cofactor vestigial-like protein 2 Short name= Vgl-2 Alternative name(s): Protein VITO1
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.