



# Recombinant Human Lethal (3)malignant brain tumor-like protein 4 (L3MBTL4)

<b>Product Code</b>	CSB-BP836703HU
<b>Abbreviation</b>	L3MBTL4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NA19
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKQPNRKRKL NMDSKERLDQ DGRLEQAEEE KPKDSTTPL SHVPSAAAQG AWSWEWYLKE QKAVAAPVEL FSKDQSFPEH ENGFQIGMRL EGIDPRHPSV FCVLSVAEVC GYRLRLHFDG YLSCYDFWTN AGSPDIHPVG WCEKTKHELH IPKGYRKDKF VWMDYLKACK LQNAPKKLFR NRSPNGPMSK EFQVGMKLEA VDRKNPSLVC VATIADIVED RLLVHFDNWD DSYDYWCDVN SPYVQPVGWC QENGRTLIAQ QGYPNPENFS WTEYLEATQT NAVPAKVFKM RLPHGFLPNM KLEVVDKRNP RLIRVATIVD VDDQRVKVHF DGWDHKYDYW VEADSPDIHP IGWCDVTGHP LEVPQRTNDL KILPGQAVCP TPGCRGIGHI RGPYSGHHS AFGCPYSMDN LKKEATLHDR LREQTQANLE SDSSHKS KSKS LCSLNFNGKH EKVNSQPRLV QQAKCLKIKG KEDIDLNL F RVLVLHPRGL EYSVEQAQQV LHQSVMSTV SAHPFRDLPL GREQHCKLLP GVADIRASQV ARWTVDEVAE FVQSLGCEE HAKCFKKEQI DGKAFLLLTQ TDIVKVMKIK LGPALKIYNS ILMFRHSQEL PEEDIASGQE VRG
<b>Source</b>	Baculovirus
<b>Target Names</b>	L3MBTL4
<b>Protein Names</b>	Recommended name: Lethal(3)malignant brain tumor-like protein 4 Short name= H-l(3)mbt-like protein 4 Short name= L(3)mbt-like protein 4
<b>Expression Region</b>	1-623
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.