



Recombinant Human Lethal (3)malignant brain tumor-like protein 4 (L3MBTL4)

Product Code	CSB-MP836703HU
Abbreviation	L3MBTL4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NA19
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKQPNRKRKL NMDSKERLDQ DGRLEQAEEE KPKDSTTPL SHVPSAAAQG AWSWEWYLKE QKAVAAPVEL FSKDQSFPEH ENGFQIGMRL EGIDPRHPSV FCVLSVAEVC GYRLRLHFDG YLSCYDFWTN AGSPDIHPVG WCEKTKHELH IPKGYRKDKF VWMDYLKACK LQNAPKKLFR NRSPNGPMSK EFQVGMKLEA VDRKNPSLVC VATIADIVED RLLVHFDNWD DSYDYWCDVN SPYVQPVGWC QENGRTLIAQ QGYPNPENFS WTEYLEATQT NAVPAKVFKM RLPHGFLPNM KLEVVDKRNP RLIRVATIVD VDDQRVKVHF DGWDHKYDYW VEADSPDIHP IGWCDVTGHP LEVPQRTNDL KILPGQAVCP TPGCRGIGHI RGPYRSGHHS AFGCPYSMDN LKKEATLHDR LREQTQANLE SDSSHKS KSKS LCSLNFNGKH EKVNSQPRLV QQAKCLKIKG KEDIDLNLF RVLVLHPRGL EYSVEQAQQV LHQSVMSTV SAHPFRDLPL GREQHCKLLP GVADIRASQV ARWTVDEVAE FVQSLGCEE HAKCFKKEQI DGKAFLLLTQ TDIVKVMKIK LGPALKIYNS ILMFRHSQEL PEEDIASGQE VRG
Source	Mammalian cell
Target Names	L3MBTL4
Protein Names	Recommended name: Lethal(3)malignant brain tumor-like protein 4 Short name= H-l(3)mbt-like protein 4 Short name= L(3)mbt-like protein 4
Expression Region	1-623
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.