



Recombinant Human Coiled-coil domain-containing protein 63 (CCDC63)

Product Code	CSB-EP836704HU
Abbreviation	CCDC63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NA47
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSVLKKNRRK DSDTPQEPSE KAKEQQAEAE LRKLRQQFRK MVESRKSFKF RNQQKIASQY KEIKTLKTEQ DEITLLL SLM KSSRNMRSE KNYMELRLLL QTKEDYEALI KSLKVLLAEL DEKILQMEKK IANQKQIFAK MQEANNPRKL KQQIHILETR LNLVTVHFDK MLTTNAKLRK EIEDLRFEKA AYDNVYQQQLQ HCLLMEKKTMLAIEQSSQA YEQRVEAMAR MAAMKDRQKK DTSQYNLEIR ELERLYAHES KLKSFLLVKL NDRNEFEEQA KREEALKAKK HVKKNRGESF ESYVAHLRL LKLAESGNLN QLIEDFLAKE EKNFARFTYV TELNNDMEMM HKRTQRIQDE IILLRSQQKL SHDDNHSVLR QLEDKLRKTT EADMYESKY GEVSKTLDLL KNSVEKLFKK INCDATKILV QLGETGKVTD INLPQYFAII EKKTNDLLLL ETYRRILEVE GAEAEIPPPF INPFWGGSSAL LKPPEPIKVI PPVLGADPFS DRLDDVEQPL DHSSLRQLVL DNYILKENRS KEVRGDSLPE KVDDFRSRKK VTM
Source	E.coli
Target Names	CCDC63
Protein Names	Recommended name: Coiled-coil domain-containing protein 63
Expression Region	1-563
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.