



# Recombinant Human Coiled-coil domain-containing protein 63 (CCDC63)

<b>Product Code</b>	CSB-YP836704HU
<b>Abbreviation</b>	CCDC63
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NA47
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSVLKKNRRK DSDTPQEPSE KAKEQQAEAE LRKLRQQFRK MVESRKSFKF RNQQKIASQY KEIKTLKTEQ DEITLLL SLM KSSRNMRSE KNYMELRLLL QTKEDYEALI KSLKVLLAEL DEKILQMEKK IANQKQIFAK MQEANNPRKL KQQIHILETR LNLVTVHFDK MLTTNAKLRK EIEDLRF EKA AYDNVYQQQLQ HCLLMEKKTMLAIEQSSQA YEQRVEAMAR MAAMKDRQKK DTSQYNLEIR ELERLYAHES KLKSFLLVKL NDRNEFEEQA KREEALKAKK HVKKNRGESF ESYVAHLRL LKLAESGNLN QLIEDFLAKE EKNFARFTYV TELNNDMEMM HKRTQRIQDE IILLRSQQKL SHDDNHSVLR QLEDKLRKTT EADMYESKY GEVSKTLDLL KNSVEKLFKK INCDATKILV QLGETGKVTD INLPQYFAII EKKTNDLLLL ETYRRILEVE GAEAEIPPPF INPFWGGG SAL LKPPEPIKVI PPVLGADPFS DRLDDVEQPL DHSSLRQLVL DNYILKENRS KEVRGDSLPE KVDDFRSRKK VTM
<b>Source</b>	Yeast
<b>Target Names</b>	CCDC63
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 63
<b>Expression Region</b>	1-563
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.