



# Recombinant Human Protein FAM47B (FAM47B)

<b>Product Code</b>	CSB-EP836705HU-B
<b>Abbreviation</b>	FAM47B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NA70
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGDRRPQDRP RSQGMDSKPW YCDKPPSKYF AKRKHRRLLRF PPVDTQNWVF VTEGMDDFRY ACQSPEDTLV CRRDEFLLPK ISLRGPQADR KSRKKKLLKK AALFSELSPV QPARKAFVEE VEAQLMTKHP LAMYPNLGKD MPPDLLLQVL KQLDPERKLE DAWARCEARE KTTEVPTESG KYPCGESCPR PPETPVSRRLR PQLPKTPVSS RRPEPPKTRV SSLRPEPPKT RVSSLHPEPP ETRASHLRVD PPETGVSHLC PEPPKTLVSS VHPEPPDTGA SHLCPEPPET RVSHLHPEPP ETGVSHLRPE PSKTQVSSLC PEPPEAGVSH LCLEPPNTHR VSSFLLQVLK LDSEKKLEDA RARCEGQEMT TEELTKPGKY HFWESCP RPF ESRMPHLRLV LPITRRMASL CLKPPKTRRV SSLCPEPTKT GASHLKELFQ EDTPSTMECV SDSLQRRHTS RKL RDFKWAG DLGVNEESIS SLFDFTPECR TTDQDQKIKK ANECASRLMY GMELDDMDEV EFLRIKYWDR RRAAPHYSYS AQRGRIRYGP WYFEPKLGKK LRSDEPLIDP KPVLEKPDEP DILDGLYGPI AFKDFILSKG YRMPGVIEKL FAKKGWTYDS VKTPIQRAVQ VYKYKEDVTD ASKED
<b>Source</b>	E.coli
<b>Target Names</b>	FAM47B
<b>Protein Names</b>	Recommended name: Protein FAM47B
<b>Expression Region</b>	1-645
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.