



Recombinant Human Protein FAM47B (FAM47B)

Product Code	CSB-EP836705HU
Abbreviation	FAM47B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NA70
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGDRRPQDRP RSQGMDSKPW YCDKPPSKYF AKRKHRRLLRF PPVDTQNWVF VTEGMDDFRY ACQSPEDTLV CRRDEFLLPK ISLRGPQADR KSRKKKLLKK AALFSELSPV QPARKAFVEE VEAQLMTKHP LAMYPNLGKD MPPDLLLQVL KQLDPERKLE DAWARCEARE KTTEVPTESG KYPCGESCPR PPETPVSRLR PQLPKTPVSS RRPEPPKTRV SSLRPEPPKT RVSSLHPEPP ETRASHLRVD PPETGVSHLC PEPPKTLVSS VHPEPPDTGA SHLCPEPPET RVSHLHPEPP ETGVSHLRPE PSKTQVSSLC PEPPEAGVSH LCLEPPNTHR VSSFLLQVLK LDSEKKLEDA RARCEGQEMT TEELTKPGKY HFWESCP RPF ESRMPHLRLV LPITRRMASL CLKPPKTRRV SSLCPEPTKT GASHLKELFQ EDTPSTMECV SDSLQRRHTS RKL RDFKWAG DLGVNEESIS SLFDFTPECR TTDQDQKIKK ANECASRLMY GMELDDMDEV EFLRIKYWDR RRAAPHYSYS AQRGRIRYGP WYFEPKLGKK LRSDEPLIDP KPVLEKPDEP DILDGLYGPI AFKDFILSKG YRMPGVIEKL FAKKGWTYDS VKTPIQRAVQ VYKYKEDVTD ASKED
Source	E.coli
Target Names	FAM47B
Protein Names	Recommended name: Protein FAM47B
Expression Region	1-645
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.