



# Recombinant Human Nitric oxide-associated protein 1 (NOA1)

<b>Product Code</b>	CSB-YP836715HU
<b>Abbreviation</b>	NOA1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NC60
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLPARLPFRL LSLFLRGSAP TAARHGLREP LLERRCAAAS SFQHSSSLGR ELPYDPVDTE GFGEGGDMQE RFLFPEYILD PEPQPTREKQ LQELQQQQEE EERQRQQRRE ERRQQNLRAR SREHPVVGHP DPALPPSGVN CSGCGAELHC QDAGVPGYLP REKFLRTAEA DGGLARTVCQ RCWLLSHHRR ALRLQVSREQ YLELVSAALR RPGPSLVLYM VDLLDLPDAL LPDLPALVGP KQLIVLGNKV DLLPQDAPGY RQRLRERLWE DCARAGLLLA PGHQGPQRPV KDEPQDGENP NPPNWSRTVV RDVRLISAKT GYGVEELISA LQRSWRYRGD VYLVGATNAG KSTLFNTLLE SDYCTAKGSE AIDRATISPW PGTTLNLLKF PICNPTPYRM FKRHRQLKGD STQAEEDLSE QEQNQLNVLK KHGYVVGRVG RTFLYSEEQK DNIPFEFDAD SLAFDMENDP VMGTHKSTKQ VELTAQDVKD AHWFYDTPGI TKENCILNLL TEKEVNIVLP TQSIVPRTFV LKPGMVLFLG AIGRIDFLQG NQSAWFTVVA SNILPVHITS LDRADALYQK HAGHTLLQIP MGGKERMAGF PPLVAEDIML KEGLGASEAV ADIKFSSAGW VSVTPNFKDR LHLRGYTPEG TVLTVRPPLL PYIVNIKGQR IKKSVAYKTK KPPSLMYNVR KKKGKINV
<b>Source</b>	Yeast
<b>Target Names</b>	NOA1
<b>Protein Names</b>	Recommended name: Nitric oxide-associated protein 1
<b>Expression Region</b>	1-698
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.