



Recombinant Human F-box/LRR-repeat protein 13 (FBXL13)

| | |
|--------------------------|---|
| Product Code | CSB-BP836729HU |
| Abbreviation | FBXL13 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8NEE6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MTPELMIKAC SFYTGHLVKT HFCTWRDIAR TNENVVLAEK MNRAVTCYNF RLQKSVFHHW HSYMEDQKEK LKNILLRIQQ IYCHKLTII LTKWRNTARH KSKKKEDELI LKHELQLKKW KNRLILKRAA AEESNFPERS SSEVFLVDET LKCDISLLPE RAILQIFFYL SLKDVICGQ VNHAWMLMTQ LNSLWNAIDF SSVKNVIPDK YIVSTLQRWR LNVLRLNFRG CLLRPKTFRS VSHCRNLQEL NVSDCPTFTD ESMRHISEGC PGVLCLNLSN TTITNRTMRL LPRHFHNLQN LSLAYCRRFT DKGLQYLNLG NGCHKLIYLD LSGCTQISVQ GFRYIANSCT GIMHLTINDM PTLTDNCVKA LVEKCSRITS LVFTGAPHIS DCTFRALSAC KLRKIRFEGN KRVTDA SFKF IDKNYPNL SH IYMADCKGIT DSSLRSLSPL KQLTVLNLAN CVRIGDMGLK QFLDGPASMR IRELNLSNCV RLSDASVMKL SERCPNLNYL SLRNCEHLTA QGIGYIVNIF SLVSI DLSGT DISNEGLNVL SRHKKLKELS VSECYRITDD GIQAFCKSSL ILEHLDVSYC SQLSDMIIKA LAIYCINLTS LSIAGCPKIT DSAMEMLSAK CHYLHILDIS GCVLLTDQIL EDLQIGCKQL RILKMQYCTN ISKKAQRMS SKVQQQEYNT NDP PRWFGYD REGNPVTELD NITSSKGALE LTVKKSTYSS EDQAA |
| Source | Baculovirus |
| Target Names | FBXL13 |
| Protein Names | Recommended name: F-box/LRR-repeat protein 13 Alternative name(s): F-box and leucine-rich repeat protein 13 |
| Expression Region | 1-735 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.