



Recombinant Human F-box/LRR-repeat protein 13 (FBXL13)

Product Code	CSB-EP836729HU-B
Abbreviation	FBXL13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEE6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTPELMIKAC SFYTGHLVKT HFCTWRDIAR TNENVVLAEK MNRAVTCYNF RLQKSVFHHW HSYMEDQKEK LKNILLRIQQ IYCHKLTII LTKWRNTARH KSKKKEDELI LKHELQLKKW KNRLILKRAA AEESNFPERS SSEVFLVDET LKCDISLLPE RAILQIFFYL SLKDVICGQ VNHAWMLMTQ LNSLWNAIDF SSVKNVIPDK YIVSTLQRWR LNVLRLNFRG CLLRPKTFRS VSHCRNLQEL NVSDCPTFTD ESMRHISEGC PGVLCLNLSN TTITNRTMRL LPRHFHNLQN LSLAYCRRFT DKGLQYLNLG NGCHKLIYLD LSGCTQISVQ GFRYIANSCT GIMHLTINDM PTLTDNCVKA LVEKCSRITS LVFTGAPHIS DCTFRALSAC KLRKIRFEGN KRVTDA SFKF IDKNYPNL SH IYMADCKGIT DSSLRSLSPL KQLTVLNLAN CVRIGDMGLK QFLDGPASMR IRELNLSNCV RLSDASVMKL SERCPNLNYL SLRNCEHLTA QGIGYIVNIF SLVSI DLSGT DISNEGLNVL SRHKKLKELS VSECYRITDD GIQAFCKSSL ILEHLDVSYC SQLSDMIIKA LAIYCINLTS LSIAGCPKIT DSAMEMLSAK CHYLHILDIS GCVLLTDQIL EDLQIGCKQL RILKMQYCTN ISKKAQRMS SKVQQQEYNT NDP PRWFGYD REGNPVTELD NITSSKGALE LTVKKSTYSS EDQAA
Source	E.coli
Target Names	FBXL13
Protein Names	Recommended name: F-box/LRR-repeat protein 13 Alternative name(s): F-box and leucine-rich repeat protein 13
Expression Region	1-735
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.