



Recombinant *Corynebacterium glutamicum* Lipoyl synthase (lipA)

Product Code	CSB-YP836806DWY
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NNJ0
Product Type	Recombinant Protein
Immunogen Species	<i>Corynebacterium glutamicum</i> (strain ATCC 13032 / DSM 20300 / JCM 1318 / LMG 3730 / NCIMB 10025)
Purity	>85% (SDS-PAGE)
Sequence	MTIAPEGRRL LRVEARNSET PIETKPRWIR NQVKNGPEYQ DMKERVAGAS LHTVCQEAGC PNIHECWESR EATFLIGGAN CSRRCDFMI NSARPEPLDR GEPLRVAESV REMQLNYSTI TGVTRDDLDD EGAWLYSEVV RKIHELNPHT GVENLVPDFS GKKDLLQEVF ESRPEVFAHN VETVPRIFKR IRPAFRYERS LDVIRQARDF GLVTKSNLIL GMGETKEEIT EALQDLHDAG CDIITITQYL RPGPLFHPIE RWVKPEEFLE HADAAKEMGF AAVMSGPLVR SSYRAGRLYA QAMEFRGEEI PAHLAHLKDT SGGSTAQEAS TLLERYGASE DTPVVVSN
Source	Yeast
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.