



Recombinant *Xanthomonas campestris* pv. *campestris* Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-BP836933XAY
Abbreviation	psd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8P7Q6
Product Type	Recombinant Protein
Immunogen Species	<i>Xanthomonas campestris</i> pv. <i>campestris</i> (strain ATCC 33913 / DSM 3586 / NCPPB 528 / LMG 568 / P 25)
Purity	>85% (SDS-PAGE)
Sequence	MSLVTSPTYV LPHRLLSSLA RALAYSKSPA TKQWLIDTVT RKFGVDLSEA QEPDPRVYPT FNAFFTRALK PGARVPDADP QALLMPADGR ISQLGPIENG RIFQAKGQSF TAAELLGDES AAVPFHNGLF ATVYLSPKDY HRVHMPWSGT LRETVHVPGR LFSVGPDAVR NVPRLFARNE RLVCHFDTDF GPMA SVMVGA LLVSGVETVW SGVEIPRYGD RITRKDYRGK GITLERFAEM ARFNYG
Source	Baculovirus
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.