



Recombinant Mouse Vacuolar protein sorting-associated protein 37B (Vps37b)

Product Code	CSB-BP837114MO
Abbreviation	Vps37b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R0J7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MAGAVSEARF AGLSLMQLHE LLEDDAQLGD MVRGMEEAQT VQLNKEMTLA SNRSLAEGNL LYQPQLDAQK ARLTQKYQEL QVLFEAYQIK KTKLSDKQSNN ASLETLALL QAEGAKIEED TENMAEKFLD GELPLDSFID VYQSKRKLAH MRRVKVEKLQ ELVLKGQRHP QAGAPPPPRV PEPSPATALP YPSLEATGLP SVVPRRIPPP PPPVPAGHVA TPFAAAMGSG QVSAYPGLQC PPLPPRVGLP SQQGFSAQLV SPYPPALPQR PPRMAPHQF GFILQ
Source	Baculovirus
Target Names	Vps37b
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 37B Short name= Vps37B Alternative name(s): ESCRT-I complex subunit VPS37B
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.