



Recombinant Mouse Angiopoietin-related protein 6 (Angptl6)

Product Code	CSB-YP837116MO
Abbreviation	Angptl6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R0Z6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p> ARCRVT LVLSPQKATS AVCRSSEATQ DSELATLRMR LGRHEELLRA LQRRAAEGGA LADEVRALRE HSLTLNTRLG QLRAQLQQEA RAEPDLGAEP AAALGLLAER ALDAEAEARR TTARLQQLDA QLREHAQLMS QHSSLLGRLQ RACAGPERGQ QQVLPLPLAP LVPLSLVGSA SNTSRRLDQT PEHQREQSLR QQGPPSSLLP TGH LAVPTRP VGPWRDCAEA HGAGHWQSGV YDLRLGRRVV AVWCEQQQEG GGWTVIQRRQ DGSVNFFTNW QHYKAGFRP EGEYWLGLEP VHQVTSRGDH ELLILLEDWG GRAARAHYDS FSLEPESDHY RLRLGQYHGD AGDSL SWHND KPFSTVDRDR DSYSGNCALY HRGGWWYHAC AHSNNGVWY HGGHYRSRYQ DGVYWAEFRG GAYSLKKAVM LTRLVRL </p>
Source	Yeast
Target Names	Angptl6
Protein Names	Recommended name: Angiopoietin-related protein 6 Alternative name(s): Angiopoietin-like protein 6 Angiopoietin-related growth factor
Expression Region	25-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.