



# Recombinant Mouse Rab11 family-interacting protein 5 (Rab11fip5)

<b>Product Code</b>	CSB-BP837136MO
<b>Abbreviation</b>	Rab11fip5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8R361
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALVRDPEPA AGSSRWLPTH VQVTVLRASG LRGKSSGAGS TSDAYTVIQV GREKYSTSVV EKTQGCPEWC EECSFELPPG ALDGLLRAQE ADAGPAPWAS GPNAACELVL TTMHRSLIGV DKFLGRATVA LDEVFRAGRA QHTQWYRLHS KPGKKEKERG EIQVTIQFTR NNLSASMF DL SMKDKPRSPF SKLKDRVKGK KKYDLESASA ILPSSALEDP ELGSLGKMGK AKGFFLRNKL RKSSLTQSNT SLGSDSTLSS TSGSLVYQGP GAELLTRSPS HSSWLSTEGG RDSIQSPKLL THKRTYSDEA SQLRAAPPRA LLELQGHLDG ASRSSLCVNG SHVYNEEPQP PLRHRSSISG PFPSSSLHS VPPRSSEEGS RSSDDSWGRG SHGTSSSEAV PGQEELSKQA KGASCSGEEE GARLPEGKPV QVATPMVASS EAVAAEKDRK PRMGLFHFFF HQGLSRSEQG RRGVSGEKGS PSLGASPHHS STGEEKAKSS WFGLRESKEP TQKPSHPVK PLTAAPVEAS PDRKQPRTSL STALSSGLER LKTVTSGGIQ SVLPASQLGS SVDTKRPKDS AVLDQSAKYY HLTHDELIGL LLQRERELSQ RDEHVQELES YIDRLLVRIM ETSPTLLQIS PGPPK
<b>Source</b>	Baculovirus
<b>Target Names</b>	Rab11fip5
<b>Protein Names</b>	Recommended name: Rab11 family-interacting protein 5 Short name= Rab11-FIP5 Alternative name(s): Rab11-interacting protein Rip11
<b>Expression Region</b>	1-645
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.