



Recombinant Mouse Rab11 family-interacting protein 5 (Rab11fip5)

Product Code	CSB-MP837136MO
Abbreviation	Rab11fip5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R361
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MALVRDPEPA AGSSRWLPTH VQVTVLRASG LRGKSSGAGS TSDAYTVIQV GREKYSTSVV EKTQGCPEWC EECSFELPPG ALDGLLRAQE ADAGPAPWAS GPNAACELVL TTMHRSLIGV DKFLGRATVA LDEVFRAGRA QHTQWYRLHS KPGKKEKERG EIQVTIQFTR NNLSASMF DL SMKDKPRSPF SKLKDRVKGK KKYDLESASA ILPSSALEDP ELGSLGKMGK AKGFFLRNKL RKSSLTQSNT SLGSDSTLSS TSGSLVYQGP GAELLTRSPS HSSWLSTEGG RDSIQSPKLL THKRTYSDEA SQLRAAPPRA LLELQGHLDG ASRSSLCVNG SHVYNEEPQP PLRHRSSISG PFPSSSLHS VPPRSSEEGS RSSDDSWGRG SHGTSSEAV PGQEELSKQA KGASCSGEEE GARLPEGKPV QVATPMVASS EAVAAEKDRK PRMGLFHFFF HQGLSRSEQG RRGSVGEKGS PSLGASPHHS STGEEKAKSS WFG LRESKEP TQKPSHPVK PLTAAPVEAS PDRKQPRTSL STALSSGLER LKTVTSGGIQ SVLPASQLGS SVDTKRPKDS AVLDQSAKYY HLTHDELIGL LLQRERELSQ RDEHVQELES YIDRLLVRIM ETSPTLLQIS PGPPK
Source	Mammalian cell
Target Names	Rab11fip5
Protein Names	Recommended name: Rab11 family-interacting protein 5 Short name= Rab11-FIP5 Alternative name(s): Rab11-interacting protein Rip11
Expression Region	1-645
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.