



Recombinant *Oryza sativa* subsp. *japonica* Probable U3 small nucleolar RNA-associated protein 11 (Os01g0810000, LOC_Os01g59500)

Product Code	CSB-BP837300OFG
Abbreviation	Os01g0810000, LOC_Os01g59500
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8S1Z1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSSLRNAIQR RAHKERAQPE SRKKFGLLEK HKDYIVRAKA FHQKEETIRK LKEKASFRNP DEFYFKMINS KTVDGIHRPK PEANKYTEDE LMLLKTDMG YILQGIQSEK KKIEKLSSML HELDNKRPNK HVYFAEDREE VKEIQSRIEQ KSSSLGLDNI PSRIKRKTAS SYRELEERKQ RVQKLEKLYA DMALQKELKK PGRKRKLRED EIENQTSRPV YKWRAQRKR
Source	Baculovirus
Target Names	Os01g0810000
Protein Names	Recommended name: Probable U3 small nucleolar RNA-associated protein 11 Short name= U3 snoRNA-associated protein 11
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.