



Recombinant Arabidopsis thaliana Putative F-box protein At2g16220 (At2g16220)

Product Code	CSB-EP837306DOA
Abbreviation	At2g16220
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8S8C7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MNSHFLTNDL ILEVLSRLPL KSVARFHCVS KRWASMFVGGSP YFKELFLTRS SAKPRLLFAI VQNGVWRFFS SPRLEKSSST LVATAEFHMK LSPNNLRIYH DNTPRYFSIG YASGLIYLYG DRYEATPLIC NPNTGRYTL PKCYTYRKAF SFFGFDPIDK QYKALSMIYP SGPGHKILT FGDGDMNWKK IKYRVLHDIY SQGICINGVL YYLGDTSDWD NDHDVTSGNV LVCFDLRSES FTFIGLECGQ LINYKGKLAV IFWDDVDDDD VKDDAIDEMH VVWLEDVEKK EWSKYAYTWT EDKLYRSQVS VVGMTASGEI VFMRKYTSE QPFYVYFNP ERNNLRRVEI QGFEAFTKFG TVYTFVDHVE DLNVYNLKHV KSVHPPSVQP EYDESDSESE EDREIII
Source	E.coli
Target Names	At2g16220
Protein Names	Recommended name: Putative F-box protein At2g16220
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.