



Recombinant Encephalitozoon cuniculi Probable serine/threonine-protein kinase RIO1 homolog (RIO1)

| | |
|--------------------------|--|
| Product Code | CSB-YP837384EKH |
| Abbreviation | RIO1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8SVI7 |
| Product Type | Recombinant Protein |
| Immunogen Species | Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKETRAEKRRK KDKSDRATVD KVLDRKRTLKV LERLQARGKL VNLQGSLCTG KESNVYLGEA STSLCSKFIK NRYSVTEEPG REGQIVPVVV KIFKTSIMSF RDRERYIRSE KRFQRFCTSN SRKLIKVWAE KEVRNLKRLN NAGIPSEPEI YLKNNILVMT QIGRCSEVAP RLRDASIKDL EGCYQQCVKI IRDMYKKAGL VHADLSEFNL LYFEGVVYVI DVGQSVEIDH DNAQRFLIMD INNINSFFSR KGVSVAKGND LFEEISGNVI PLYLKDIDIG RDAFIPSRVS EVGNEEDLLA FAADSRREF GSTTDSLSS TGEASVEDSD ASLGATESRG GGNIREEKRR KKKTVKELNR IRRASRISKK EKKRIFKRYI GMKKRKN |
| Source | Yeast |
| Target Names | RIO1 |
| Protein Names | Recommended name: Probable serine/threonine-protein kinase RIO1 homolog EC= 2.7.11.1 |
| Expression Region | 1-387 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.