



Recombinant *Drosophila melanogaster* COP9 signalosome complex subunit 3 (CSN3)

Product Code	CSB-MP837396DLU
Abbreviation	CSN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SYG2
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGSALENYVN QVRTLSASGS YRELAELPE SLSLLARNWS ILDNVLETLD MQQHSLGVLY VLLAKLHSAS TANPEPVQLI QLMRDFVQRN NNEQLRYAVC AFYETCHLFT EFVVQKNLSI LGIRIISRAI DQIRQLETQL TPIHADLCLL SLKAKNFSVV LPYLDADITD ISTVAAECKT QQQQQSQHAD ANNDAKYFLL YFYYGGMIYT AVKNYERALLY FFEVCITTPA MAMSHIMLEA YKKFLMVSLI VEGKIAYIPK NTQVIGRFMK PMANHYHDLV NVYANSSSEE LRIILKYSE AFTRDNNMGL AKQVATSLYK RNIQRLTKTF LTLSDVAS RVQLASAVEA ERYILNMIKS GEIYASINQK DGMVLFKDDP EKYNPEMFL NVQNNITHVL DQVRQINKME EEIILNPMYV KKALGSQDDD LTSQHPTFS GDPTD
Source	Mammalian cell
Target Names	CSN3
Protein Names	Recommended name: COP9 signalosome complex subunit 3 Short name= Dch3 Short name= Signalosome subunit 3
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.