



Recombinant Human Guanine nucleotide exchange factor for Rab-3A (RAB3IL1)

Product Code	CSB-YP837424HU
Abbreviation	RAB3IL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TBN0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MWSGPPQPDQ GLPPPLAAVP VPWKSTDPCQ GHRESPGALV ETSAGEEAQG QEGPAAAQLD VLRLRSSME IREKGSEFLK EELHRAQKEL KLKDEECERL SKVREQLEQE LEELTASLFE EAHKMVREAN MKQAASEKQL KEARGKIDML QAEVTALKTL VITSTPASPEN RELHPQLLSP TKAGPRKGHS RHKSTSSTLC PAVCPAAGHT LTPDREGKEV DTILFAEFQA WRESPTLDKT CPFLERVYRE DVGPCLDFTM QELSVLVRAA VEDNTLTIEP VASQTLPTVK VAEVDCSSTN TCALSGLTRT CRHRIRLGDS KSHYYISPSS RARITAVCNF FTYIRYIQQG LVRQDAEPMF WEIMRLRKEM SLAKLGFFPQ EA
Source	Yeast
Target Names	RAB3IL1
Protein Names	Recommended name: Guanine nucleotide exchange factor for Rab-3A Alternative name(s): Rab-3A-interacting-like protein 1 Short name= Rab3A-interacting-like protein 1 Rabin3-like 1
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.