



# Recombinant Human Protein phosphatase Slingshot homolog 3 (SSH3)

<b>Product Code</b>	CSB-EP837447HU
<b>Abbreviation</b>	SSH3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TE77
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ALVTVSRSP PGSGASTPVG PWDQAVQRRS RLQRRQSFVAV LRGAVLGLQD GGDNDDAAEA SSEPTEKAPS EEELHGDQTD FGQGSQSPQK QEEQRQHLHL MVQLLRPQDD IRLAAQLEAP RPPRLRYLLV VSTREGEGLS QDETVLLGVD FPDSSSPSCT LGLVLPLWSD TQVYLDGDGG FSVTSGGQSR IFKPISIQTM WATLQVLHQA CEALGSGLV PGGSSALTWAS HYQERLNSEQ SCLNEWTAMA DLESLRPPSA EPGGSSESEQ MEQAIRAELW KVLVDVSDLES VTSKEIRQAL ELRLGLPLQQ YRDFIDNQML LLVAQRDRAS RIFPHLYLGS EWNAAANLEEL QRNRVTHILN MAREIDNFYP ERFTYHNVRL WDEESAQLLP HWKETHRFIE AARAQGTHVL VHCKMGVSRS AATVLAYAMK QYECSSLEQAL RHVQELRPIA RPNPGFLRQL QIYQGILTAS RQSHVWEQKV GGVSPPEEHPA PEVSTPFPPPL PPEPEGGGEE KVVGMEESSQA APKEEPGPRP RINLRGVMRS ISLLEPSLEL ESTSETSDMP EVFSSHESSH EEPLQFPFQL ARTKGGQQVD RGPQPALKSR QSVVTLQGSA VVANRTQAFQ EQEQGGQGGQ GEPGISSTPR FRKVVRQASV HDSGEEGEA
<b>Source</b>	E.coli
<b>Target Names</b>	SSH3
<b>Protein Names</b>	Recommended name: Protein phosphatase Slingshot homolog 3 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): SSH-like protein 3 Short name= SSH-3L Short name= hSSH-3L
<b>Expression Region</b>	2-659
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.