



Recombinant Human Protein phosphatase Slingshot homolog 3 (SSH3)

Product Code	CSB-YP837447HU
Abbreviation	SSH3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TE77
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ALVTVSRSP PGSGASTPVG PWDQAVQRRS RLQRRQSFVAV LRGAVLGLQD GGDNDDAAEA SSEPTEKAPS EEELHGDQTD FGQGSQSPQK QEEQRQHLHL MVQLLRPQDD IRLAAQLEAP RPPRLRYLLV VSTREGEGLS QDETVLLGVD FPDSSSPSCT LGLVLPLWSD TQVYLDGDGG FSVTSGGQSR IFKPISIQTM WATLQVLHQA CEALGSGLV PGGSSALTWAS HYQERLNSEQ SCLNEWTAMA DLESLRPPSA EPGGSSESEQ MEQAIRAEW KVLVDVSDLES VTSKEIRQAL ELRLGLPLQQ YRDFIDNQLL LLVAQRDRAS RIFPHLYLGS EWNANLEEL QRNRVTHILN MAREIDNFYP ERFTYHNVRL WDEESAQLLP HWKETHRFIE AARAQGTHVL VHCKMGVSRS AATVLAYAMK QYECSSLEQAL RHVQELRPIA RPNPGFLRQL QIYQGILTAS RQSHVWEQKV GGVSPPEHPA PEVSTPFPPPL PPEPEGGGEE KVVGMEESSQA APKEEPGPRP RINLRGVMRS ISLLEPSLEL ESTSETSDMP EVFSSHESSH EEPLQFPFQL ARTKGGQQVD RGPQPALKSR QSVVTLQGSA VVANRTQAFQ EQEQGGGGQQ GEPGISSTPR FRKVVRQASV HDSGEEGEEA
Source	Yeast
Target Names	SSH3
Protein Names	Recommended name: Protein phosphatase Slingshot homolog 3 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): SSH-like protein 3 Short name= SSH-3L Short name= hSSH-3L
Expression Region	2-659
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.