



# Recombinant Pyrococcus furiosus Peptide chain release factor subunit 1 (prf1)

<b>Product Code</b>	CSB-MP837591FHW
<b>Abbreviation</b>	prf1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8U0J4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pyrococcus furiosus (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRVMTSRDAQ LYELKKKIDE LKKIRGRGTE LISLYIPAGY DLSKVMQQLR EEYSTAQNIK SKTTRKNVLG ALERAMQHLK LYKQTPENGL ALFVGNVSEM EGVTDIRLWA IIPPEPLNVR LYRCDQTFVT EPLEEMLRVK EAYGLITVEK NEATIGILRG KRIEVLDEL TSNVPGKTRAG GQSARRYERI REQETHEFMK RIGEHRANRIF LPLLESGELK GIIVGGPGPT KEEFVEGDYL HHELKKKIIG IVDISYHGEY GLRELVAKAA DILRDHEVIR ERNLVNEFLK HIVKDTGLAT YGEREVRKAL ELGAVDTLLI SEG YDKVRVH VKCNCGWEE LKTMSEEEYE AYKKRIQTCP KCGSQNLTFE KWEVAEELIK IAEEAGSNVE IISLDTEEGQ QFYRAFGGLG AILRFRIQGV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	prf1
<b>Protein Names</b>	Recommended name: Peptide chain release factor subunit 1 Alternative name(s): Translation termination factor aRF1
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.