



Recombinant Human ETS domain-containing protein Elk-4 (ELK4)

Product Code	CSB-MP007607HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28324
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDSAITLWQF LLQLLQKPQN KHMICWTSND GQFKLLQAEV VARLWGIRKN KPNMNYDKLS RALRYYYVKN IIKKVNGQKF VYKFVSYPEI LNMDPMTVGR IEGDCESLNF SEVSSSSKDV ENGGKDKPPQ PGAKTSSRND YIHSGLYSSF TLNSLNSSNV KLFKLIKTEP PAEKLAEEKS PQEPTPSVIK FVTTPSKKPP VEPVAATISI GPSISPSSEE TIQALETLSV PKLPSLEAPT SASNVMTAFA TTPPISSIPP LQEPPTPSP PLSSHPDIDT DIDSVASQPM ELPENLSLEP KDQDSVLEK DKVNNSSRSK KPKGLELAPT LVITSSDPSP LGILSPSLPT ASLTPAFFSQ TPIILTPSPL LSSIHFWSTL SPVAPLSPAR LQGANTLRFQF PSVLNSHGPF TSLGLDGPST PGPFPDLQK T
Source	Mammalian cell
Target Names	ELK4
Protein Names	Recommended name: ETS domain-containing protein Elk-4 Alternative name(s): Serum response factor accessory protein 1 Short name= SAP-1 Short name= SRF accessory protein 1
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene is a member of the Ets family of transcription factors and of the ternary complex factor (TCF) subfamily. Proteins of the TCF subfamily form a ternary complex by binding to the the serum response factor and the serum response element in the promoter of the c-fos proto-oncogene. This protein is phosphorylated by the kinases, MAPK1 and MAPK8. Several transcript variants have been described for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.